



X-Ray Grid



Radiography System



Software



Computed Radiography



Processor




X-Ray Accessories




RF Treatment





- · Crystal clear image quality
- · Compact, lightweight and wall-mountable to fit anywhere
- Superior scanning speed and process



Dimensions

D 803mm (31.6")

W 456mm (17.9")

H 134mm (5.4")







Cassette Sizes

35x43cm, 24x30cm and 18x24cm

(default option consists of two 35x43cm)

Auto-alignment after transportation


ExamVue CR PRO | Compact & Flexible CR Solution .

Specifications

70 large (35x43cm) plates per hour Available Speeds 50 large (35x43cm) plates per hour Resolution 16-bits per pixel, 65,535 graytones - FDA510(k), CE, UL, cUL, FCC Regulatory Approvals Not approved for mammography Simple and intuitive ExamVue Acquire Software DICOM store/send to any DICOM destination Annotations and study notes and measuring tools Cost savings by eliminating film and chemistry Increase productivity Features High quality images Wall mount option for limited space Automatic cassette and plate handling Temperature: 15~30°C (59~86°F) Temperature gradient 0.5°C/min **Operating Conditions** ■ Humidity: 15~95 RH Magnetic field: max 1,260 micro tesla

Quantor⁺ Imaging Software

for Medical & Veterinary



Next-generation touch-screen user interface

Streamlined intuitive workflow

Automation using DICOM work list server

Position auto-matching

Fast, high-quality image processing designed for CR

Full image manipulation functionality

Fast, easy image verification

Simple data management

Multi-tasking capabilities for maximum efficiency

Secure, high-speed local database for facilities without an information system

Full range of image distribution & output

Support various screen resolutions for notebook or tablet PCs

JPI Healthcare Co..Ltd.

- 🐔 608-ho, 28, Digital-ro 33-gil, Guro-gu, 08377 Seoul, Korea
- www.jpi-korea.com

