




X-Ray Grid



Radiography System




Software



Computed Radiography



Processor



X-Ray Accessories



RF Treatment




Diagnostics



APRON

(X-Ray Protective Apron)



Balanced Color Combination

Unique Design

Dark Beige Khaki Dark Beige Navy Blue

2 Colors

APRON Will Protect You

Ventilation Mesh



Velcro, Buckle Type

Comfort Shoulder Pad

Easy to Clean

APRON | X-Ray Protective Apron

Specifications

 2 Colors	 D x W (cm)	Size \ Type	Velcro	Buckle	Coat	Vest	Skirt
		Pb equivalent (mmPb)	0.25 / 0.5	0.25 / 0.5	0.25	0.25 / 0.5	0.25
Small	90 X 60	90 X 60	90 X 60	-	-	-	-
Medium	100 X 60	100 X 60	100 X 60	-	-	-	-
Large	110 X 60	110 X 60	110 X 60	69 X 60	57 X 138	-	-
X Large	110 X 70	110 X 70	110 X 70	-	-	-	-
XX Large	110 X 80	110 X 80	110 X 80	-	-	-	-
Weight	Minimum 2.475 kg	Minimum 2.535 kg	Minimum 5.25 kg	Minimum 2.64 kg	Minimum 3.27 kg	-	-



DENTAL APRON

Free Size (cm)



Pb equivalent (mmPb)
0.25 / 0.5

Minimum 1.59 kg



NECK APRON

Free Size (cm)



Pb equivalent (mmPb)
0.25 / 0.5

Minimum 2.575 kg

OTHER APRONS

Pb equivalent 0.5 mmPb

-  Taffeta
-  Glove
-  Thyroid
-  Gonad
-  Ovary

JPI Healthcare Co.,Ltd.
 608-ho, 28, Digital-ro 33-gil, Guro-gu, 08377 Seoul, Korea
 82-2-2108-2580 82-2-2108-1180
 www.jpi-korea.com

