



## ExamVue DR Software for Veterinary



Since 1980



X-ray Grids



Radiography System



Detector



Software



# ExamVue DR



Computed Radiography



Integration of various flat panel detector by different manufactures



Multiple language support



### Animal Information

Modality Work List Generator provided as bundle



Fast and professional image processing for each individual examination



X-ray Accessories

### Useful Imaging and Annotation Tools

- Free Rotation, Black Blocker, Image Export, Automatic or manual removal of grid lines
- DICOM3.0 Compatible, DICOM Print, DICOM Send

### Convenient Diagnosis Function

- Freely Configurable body parts
- Fast and safe registration of emergency patients



Gastroenterology

### User Friendly Design

- Optimized GUI design for touch screen

### Simple & Fast Workflow

- Faster Image Acquisition using the preview function



## Proposed Imaging Solution for VET


### Multiple Language Support



▶ Easy process of adding unsupported languages.

### Image Optimization Features

Enhanced various Image processing technology & Parameter adjustment Easily modifiable features such as Sharpness, Edge, Contrast, Latitude, Gamma, Noise reduction Customizable image processing parameter per user's preference capability.




Hard

Medium

Soft

### Fast Image Processing Time

- Avg. image processing time of 1.3 second required for 9MPx
- Avg. image processing time of 2 seconds required for 16MPx



(Note: Workstation used for testing's specification.  
CPU: Intel Core i5-3470 (3.20GHz), RAM 8.00GB)

JPI Healthcare Co., Ltd.

608-ho, 28, Digital-ro 33-gil, Guro-gu, 08377, Seoul, Korea

+82.2-2108-2580      +82.2-2108-1180

www.jpi-korea.com

